
CARBOSTENT™

ISTHMUS

LOGIC

The main title 'ISTHMUS LOGIC' is centered on the page. 'ISTHMUS' is in a large, bold, black sans-serif font. To its right is a square graphic containing a stylized, overlapping leaf or scale pattern. Below 'ISTHMUS' is the word 'LOGIC' in a smaller, bold, black sans-serif font. Above 'ISTHMUS' is the word 'CARBOSTENT™' in a smaller, bold, black sans-serif font.

IC0711_137.004 Issued 2022/04 Rev. 7

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of **Alvimedica** group
www.alvimedica.com

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DESCRIPTION

The "ISTHMUS LOGIC CARBOSTENT" device consists of an iliac stent firmly held on the distal end of a semicompliant balloon catheter.

The peripheral iliac **stent** is a flexible implantable device that can be expanded using a PTA balloon catheter.

The stent is made of cobalt chromium alloy and is coated with rCarbofilm™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs.

Coating the substrate provides it with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Two radio-opaque platinum markers at either end of the stent allow for accurate positioning over the lesion to be treated.

The over-the-wire type **balloon catheter** provides a safe means for delivering the iliac stent to the lesion to be treated.

The body of the catheter consists of two coaxial lumens: one for balloon inflation and deflation, the other for guidewire insertion and retraction.

Two radio-opaque markers, located outside the usable length of the balloon, enable for accurate placement across the stenosis.

The proximal end of the catheter has a Y connector: the first input to the connector (marked "WIRE") gives access to the central lumen for the guidewire, while the other (marked "BALLOON") allows for balloon inflation/deflation.

The Manufacturer directly manufactures the iliac stent and carries out all quality controls both during production and on the finished product, according to the norms of Good Manufacturing Practice.

INTENDED USE

The stent is intended for use in the treatment of atherosclerotic pathologies of the peripheral vessels, particularly of the iliac arteries, in order to maintain vessel patency.

INDICATIONS

The stent is recommended for the treatment of obstructive pathologies of the peripheral vessels, particularly of the common or external iliac arteries. The stent is also indicated in the case of inadequate results (residual stenosis > 30%) or dissection of the vessel wall following PTA procedure with or without stenting.

CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations:

- Vessel to be treated with evidence of massive thrombosis prior to PTA procedure.
- PTA treatment technically unsuccessful, for example due to impossibility of accessing the site with the delivery system.
- Aneurysm immediately adjacent to the stent implant site.
- Disorders that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to pharmaceuticals and/or contrast medium commonly used in PTA procedures.
- Lesions considered untreatable with PTA or other interventional techniques.
- Lesions with severe calcification that resists dilatation with risk of vessel wall rupture.

MODEL

Each ISTHMUS LOGIC CARBOSTENT device is identified by a model code and a lot number; the available models are listed in table 1: the sales code is made up of the letters IC followed by the two letters that identify the type of device, of two or three figures that indicate the nominal diameter of expansion of the stent and by a further two figures that indicate the length of the stent followed by the letter L or S according to whether the delivery catheter is 135 cm or 75 cm long.

The lot number allows traceability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included in each box; these labels can be attached to the implant patient's medical chart.

PACKAGING

The "ISTHMUS LOGIC CARBOSTENT" device is supplied sterile and individually packed in a pouch that **must not be brought into a sterile area**.

The manufacture sterilisation uses an ethylene oxide and CO₂ mixture. Sterility is guaranteed while the package remains intact and until the expiry date printed on the packaging (USE BEFORE DATE).

WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or re-sterilise. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The "ISTHMUS LOGIC CARBOSTENT" device is designed to perform as a system. Do not use its components separately.
- The "ISTHMUS LOGIC CARBOSTENT" delivery catheter is not a peripheral dilatation catheter; it should be used for stent deployment only.
- Do not use the "ISTHMUS LOGIC CARBOSTENT" if its packaging has been opened or damaged or if the period of sterility has expired: in such cases product sterility is not guaranteed.
- The "ISTHMUS LOGIC CARBOSTENT" should be carefully handled to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce mechanical alterations.
- The ISTHMUS LOGIC CARBOSTENT should be used by physicians thoroughly trained in percutaneous transluminal angioplasty (PTA) and in peripheral and/or coronary stent placement.
- The ISTHMUS LOGIC CARBOSTENT should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.

- Use standard techniques for placement of an appropriately sized femoral or brachial sheath that is compatible with accessing vessel.
- Cross-over technique is allowed only if iliac bifurcation angle is more than 45°.
- Do not pressurise the system until the stent is bifurcated across the lesion to be treated.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially deployed stent. Attempted repositioning may cause severe vessel damage.
- Do not attempt to clean or re-sterilise devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- Follow manufacturer's indications for the use of accessories (sheath and/or guiding catheter, guidewire, haemostatic valve).
- If resistance is encountered at any time during the insertion procedure, do not use force: withdraw the system, the guidewire and the sheath as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.
- Implantation of an iliac stent may cause dissection of the vessel distal and/or proximal to implantation site and may also cause sudden occlusion of the vessel, making an additional intervention necessary (surgical vascular replacement operation, further dilatation, placement of additional stents, or other procedures).
- Extreme caution should be exercised to reduce the risk of embolism of the plaque or material of thrombotic origin.
- If the stent is lost in the vessel, recovery procedures can be initiated. The procedures can, however, cause injury to vessels and/or to the vessel access site.
- During catheter inflation do not exceed the Rated Burst Pressure.
- Never use air or any other gas to inflate the balloon.

PRECAUTIONS

- Use of a gauged inflation device is strongly recommended.
- When excessive vessel tortuosity associated with widespread atherosclerotic plaques is present, catheter advancement may be difficult. In such cases incorrect handling could cause dissection, perforation or rupture of the vessel.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 2,890 G/cm (28.90 T/m)

Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest temperature change	2.4 °C	2.9 °C

In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

DRUG REGIME

Studies in clinical literature indicate the need to administer pre-operation anti-platelet therapy, anticoagulant therapy during the procedure and post-operation anti-platelet treatment.

Given the complexity of the procedure and the many factors influencing patient condition, definition of the appropriate therapy is left to the discretion of the physician.

DIRECTIONS FOR USE

Stent + delivery catheter preparation

Prior to stent delivery, administer Heparin and check that the patient's Activated Clotting Time (ACT) is greater than 180 sec.

The stent should be implanted in the target lesion using its delivery system.

- After having examined the package for any damage, remove the "ISTHMUS LOGIC CARBOSTENT" and bring it into a sterile area.
 - Check that the delivery catheter does not present kinks, bends or other damage.
 - Pull back the distal stylet
 - Carefully slide the protective cover off the stent by grasping the cover at the distal end, then the stylet. Check that the stent is intact and that it is still to the centre of the balloon.
- Damage to the "ISTHMUS LOGIC CARBOSTENT" may impair its performance.**
- Prepare the inflation device according to the manufacturer's instructions.
 - Purge any air from the balloon the stent is mounted on as follows:
 - 1) Fill the inflation device with 4 ml of contrast mixture;
 - 2) After connecting the inflation device to the Y connector input (marked "BALLOON") of the delivery catheter, point the distal tip of the catheter downwards.
 - 3) Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to normal** as the system fills with contrast medium.
 - 4) Without introducing air, repeat step 3 and aspirate for 10-15 seconds until bubbles no longer appear.

Stent insertion

The "ISTHMUS LOGIC CARBOSTENT" devices of nominal diameter 6.0 mm are compatible with sheaths of diameter 6F (2.0 mm) or greater

The "ISTHMUS LOGIC CARBOSTENT" devices of nominal diameter 7.0, 8.0, 9.0 and 10.0 mm are compatible with sheaths of diameter 6F large lumen (2.16 mm) or greater.

Guidewires of 0.035 inches (0.89 mm) or less are compatible for use with

the "ISTHMUS LOGIC CARBOSTENT". The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- a) Keep the delivery catheter the stent is pre-mounted on at ambient pressure during insertion.
- b) Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- c) Ensure that the haemostasis valve is completely open before introducing the stent + delivery catheter system.
- d) Slowly advance the system over the guidewire until the stent reaches the treatment site.
- e) When advancing the delivery system into the target vessel, make sure that the stent and the delivery balloon move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon markers in relation to the stent markers.

If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the system, the guidewire and the sheath as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.

Stent deployment and expansion

- a) Position the stent across the target lesion.
- b) Deploy the stent by slowly inflating the balloon to the nominal diameter. The table 2 (pressure/diameter) lists the diameters of the system as balloon inflation pressure varies (column I), for the five nominal diameter groups: 6.0 mm (column II), 7.0 mm (column III), 8.0 mm (column IV), 9.0 mm (column V) and 10.0 mm (column VI). Values on dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a balloon inflation pressure of about 5 atm.

Stent elastic recoil is between 2 and 7 % * depending on stent model and expansion diameter.

*** NOTE: These values are obtained from the results of in-vitro testing**

- c) Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation.
- d) Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTA procedure.
- e) If post-dilatation is required, a PTA balloon can be used. Use extreme care:
not to expand the 6.0 mm nominal diameter stent over 9.0 mm;
not to expand the 7.0 mm nominal diameter stent over 9.0 mm;
not to expand the 8.0 mm nominal diameter stent over 9.0 mm;
not to expand the 9.0 mm nominal diameter stent over 11.0 mm;
not to expand the 10.0 mm nominal diameter stent over 11.0 mm.

The final diameter of the stent must be suitable for the reference vessel diameter.

Ensure that the stent is fully expanded.

COMPLICATIONS

Stent implantation may give rise to the following complications:

- Death
- Iliac necrosis
- Dissection, perforation, rupture of the iliac artery
- Stent embolism
- Chronic stent occlusion (restenosis)
- Vessel occlusion
- Arterial spasm
- Infection
- Access site haematoma
- Drug and contrast medium induced collateral effects (anti-coagulant/anti-platelet)
- Acute or sub-acute thrombosis localised at stent level
- Trophic alterations of the lower limbs with consequent necrosis, gangrene and possible need to amputate limb

LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, The Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. The buyer shall return the device to the Manufacturer, which reserves the right to inspect the returned device and, at its sole discretion, to determine whether the device is truly defective in manufacturing or material. The warranty consists exclusively in the replacement of the defective device with a device of the same fabrication or equivalent.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for replacement of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package.

Furthermore, the Manufacturer declines any responsibility related to the consequences

arising from the treatment choices and the methods of use or the application of the device; The Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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Pressure above Rated Burst Pressure
Do not exceed Rated Burst Pressure

These values result from in vitro testing

If postdilatation is required, a PTA balloon can be used. Use extreme care:
not to expand the 6.0 mm nominal diameter stent over 9.0 mm
not to expand the 7.0 mm nominal diameter stent over 9.0 mm
not to expand the 8.0 mm nominal diameter stent over 9.0 mm
not to expand the 9.0 mm nominal diameter stent over 11.0 mm
not to expand the 10.0 mm nominal diameter stent over 11.0 mm

Carbostent & Implantable Devices