PERIPHERAL
VASCULAR

NiTiDES

Amphilimus™ (Sirolimus + Fatty Acid)
eluting Peripheral Self-Expanding stent

abLuminal
reservoir
technology 

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of  Alvimedica Group
www.alvimedica.com

CE
0373




IC0740.137.101 Issued 202105 Rev.0

NITIDES


Amphilimus™ (Sirolimus + Fatty Acid)
eluting Peripheral Self-Expanding stent

FIGURE 1


1	Radiopaque tip
2	Guide wire lumen
3	Outer Sheath
4	Stent
5	Stent stopper
6	Coil
7	Catheter body
8	Handle
9	Thumbwheel
10	Security key (detail B: instruction for use)
11	Pull back
12	Luer hub
13	Radiopaque marker of the sheath
14	Radiopaque markers of the stent

HOW TO UNLOCK THE SYSTEM


DEPLOYMENT OF THE STENT

Turn thumbwheel

Optional quick release

English	2
Français	4
Deutsch	7
Italiano	9
Español	12
Dansk	14
Português	17
Suomi	19
Svenska	22
Nederlands	24
Ελληνικά	27
Türkçe	29
Česky	31
Eesti	34
Magyar	36
Latviski	38
Lietuvių kalba	41
Slovenčina	43
Slovenščina	45
Norsk	48
Македонски	50
Srpski	53
Hrvatski	55
Română	58
Polski	60
Български	63
Русский	65
简体中文	68
Bahasa Indonesia	70
Український	72

NITIDES

αCarbofilm™-coated self-expanding Amphillimus™ (Sirolimus+fatty acid) eluting peripheral stent system**1. DEVICE DESCRIPTION**

The NiTiDES device family is composed of a self-expanding stent and a delivery system for its correct placement and deployment.

The **stent** is a flexible, self-expanding implantable device. It is made of Nickel and Titanium alloy (Nitinol) and coated with αCarbofilm™, a thin layer of carbon with a high-density turbostratic structure substantially identical to that of the pyrolytic carbon.

Coating the substrate provides it with the biocompatible and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Six radiopaque tantalum markers are applied to the stent (three at each end) to enable correct placement at the site of the target lesion.

The outer surface of the stent has dedicated grooves, fully coated with αCarbofilm™, for containing the pharmaceutical formulation Amphillimus™, which is composed of the drug Sirolimus and a mixture of long-chain fatty acids.

The specific drug dosage for the vessel surface treated is approximately 0.9 µg/mm², which corresponds to a minimum dose of 270 µg on the smaller stent (6.0x20 mm) and a maximum dose of 2420 µg on the larger stent (8.0x150 mm).

The **delivery system** utilized is an over-the-wire type, ensuring that the stent can be positioned and released safely at the lesion. It is shown schematically in Figure 1 and consists of:

- a catheter body (7) 6F in diameter composed of:
 - a retractable outer sheath, the distal, semi-transparent section (3) which contains the stent (4); the sheath has a radiopaque collar (13) near the tip of the catheter;
 - an inner lumen comprising a coil (6) and an inner tube (2) compatible with a 0.035 inch guidewire; the inner lumen extends along the catheter from a luer fitting (12) secured proximally to the handle and terminates distally at the tip (1) in radiopaque material, so as to facilitate the location and placement of the stent when viewed on an imaging system. A stopper (5) integral with the inner lumen prevents axial displacement of the stent during retraction of the outer sheath;
- a handle (8) with a delivery system, which allows the positioning of the catheter and the deployment of the stent in situ.

The handle has an ergonomic design allows the operator to grasp it and use it even with one hand. It has a luer port (12) at the proximal end, which is used to flush the guidewire lumen, a safety key (10), a wheel (9) and a pullback tail (11). The Manufacturer manufactures the peripheral stent + deployment catheter system directly and implements all quality control procedures both during production and on the finished products in accordance with Good Manufacturing Practice guidelines.

2. INTENDED USE OF THE DEVICE

The self-expanding peripheral stent system is intended for the treatment of the Peripheral Artery Disease affecting the superficial femoral artery and the proximal popliteal artery.

3. INDICATION OF THE DEVICE

The device is indicated for the treatment of stenosis of femoropopliteal arteries above-the-knee in adult symptomatic [Rutherford category 2-5] patients. The treated lesion, of maximum length of 14 centimeters, must be either occluded or with stenosis higher than 50%, de-novo or restenotic, not previously treated with a stent.

4. CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations:

- Lesions considered untreatable with PTA or other interventional techniques.
- Women who are pregnant, breastfeeding.
- Lesions localized in the two distal thirds of the popliteal artery or at the knee joint.
- Clinical conditions that impede the use of anti-platelet and/or anticoagulant therapy.

5. WARNINGS

- The use of the device is not recommended in patients with known hypersensitivity or allergies to Sirolimus, fatty acids (such as stearic acid, palmitic acid, behenic acid) or the metal components of the stent (such as Nickel, Titanium and Tantalum).
- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The device is designed to perform as a system. Do not use its components separately.
- Do not use the device if its packaging has been opened or damaged, incorrectly stored, or if it is past the stated expiry date. In such cases product sterility is not guaranteed.
- The device should be used only by physicians trained specifically to perform percutaneous transluminal angioplasty (PTA) and peripheral stent implantation.
- If resistance is encountered at any stage of the procedure, do not force the system: withdraw the stent system and the entire delivery catheter, the guidewire and the introducer as a single unit. Applying excessive force and/or incorrect handling may result in the stent being accidentally deployed or in damage to the placement catheter.
- At any time during operation, incorrect handling could cause dissection, perforation or rupture of the vessel.

- Do not advance or withdraw the stent system during deployment: it may cause serious damage to vessel walls.
- The implantation of a self-expanding peripheral stent can cause dissection of the vessel distally and/or proximally to the stent and could also cause acute blockage of the vessel, making an additional operation necessary (surgical vascular replacement operation, further dilatation, deployment of additional stents or other procedures).
- Once the stent has been partly deployed, it cannot be retrieved using the delivery system.

6. PRECAUTIONS

- **Never use** the quick-release system at the start of stent expansion. Use it **only if necessary**, after releasing at least **half the length** of the stent.
- The device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage it or produce mechanical alterations.
- The device should not come into contact with liquids prior to preparation and positioning.
- The device's placement system is not intended for use with contrast medium injectors.
- Where stents are implanted in tandem, the relative radiopaque markers should overlap. Overlap should not involve more than two stents.
- In case of overlapping stents, less than 20 mm overlapping length is recommended.
- To avoid damaging the stent, use extreme care whenever passing a guidewire or a balloon catheter or another device through the stent that has been deployed.

7. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the manufacturer's Nitinol stent product lines are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3-Tesla (3T)
- Maximum spatial field gradient of 6,140 G/cm (61.40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of:
 - 1 W/kg (Normal Operating Mode) at 1.5T
 - 1 W/kg (Normal Operating Mode) at 3T

RF Heating

Under the scan conditions defined above, the manufacturer's Nitinol stent product lines are expected to produce a maximum temperature rise of less than 6.8 °C after 15 minutes of continuous scanning at an MR system reported, whole body averaged specific absorption rate (SAR) of 1.0 W/kg.

Single Stent Configuration

Highest temperature change

Overlapping Stent Configuration

Highest temperature change

	1.5 T	3.0 T
<i>Highest temperature change</i>	≤ 6.3 °C	≤ 6.2 °C
<i>Highest temperature change</i>	≤ 6.8 °C	≤ 6.3 °C

MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately 0.7 cm from the manufacturer's Nitinol stents when imaged with gradient-echo and spin echo pulse sequences in a 3T MRI system.

8. HOW SUPPLIED

The NiTiDES device is supplied sterile (using an ethylene oxide and CO₂ mixture), for single use, wrapped in two pouches that must not be brought into a sterile area. **NOTE: the inner pouch is the only sterile barrier.**

Sterility and chemical stability are guaranteed while the package remains intact and until the expiry date printed on the packaging (USE BEFORE DATE).

Contents: one (1) αCarbofilm™-coated self-expanding Amphillimus™ (Sirolimus+fatty acid) eluting peripheral stent system

9. STORAGE

The product should be stored as indicated on the packaging, in a cool and dry place, with storage temperature max 25 °C.

10. MODELS

The NiTiDES device is identified by a model code and a lot number; available product codes are listed in Table 1.

11. DRUG REGIME

Studies in clinical literature indicate the need to administer pre-operation anti-platelet therapy, anticoagulant therapy during the procedure and post-operation anti-platelet treatment.

Antiplatelet treatment during the period after the procedure involves the administration of:

- Thienopyridines (Ticlopidine or Clopidogrel or Prasugrel in the doses indicated in the respective instruction booklets or according to medical prescription)
 - Aspirin at a minimum amount of 100 mg/day indefinitely
- Given the complexity of the procedure and the many factors influencing patient condition, definition of the appropriate therapy is left to the discretion of the physician.

12. INTERACTIONS WITH DRUGS

Although no specific clinical data are available, certain drugs, such as Tacrolimus, which act through the same binding protein (FKBP), may interfere with the efficacy of Sirolimus. Sirolimus is metabolized by CYP3A4. Strong CYP3A4 inhibitors may cause increased exposure to Sirolimus, especially in the case of multiple implant. Systemic exposure to Sirolimus if the patient is taking concomitant systemic immunosuppressive drugs must also be considered.

13. CLINICAL TRIAL WITH NITIDES

Principal safety and effectiveness of NITIDES stent have been established from the pre-market pivotal clinical trial ILLUMINA.

Design of the study: Within a multicenter, single arm study, 100 patients with symptomatic femoropopliteal disease with a maximum lesion length of 14cm were treated with the NITIDES stent. The primary efficacy (primary patency) and primary safety endpoint (a composite of major adverse events: death, target limb amputation, target limb ischemia requiring surgical intervention or surgical repair of target vessel or clinically driven target lesion revascularization, worsening of the Rutherford score by 2 classes, or to class 5 or 6) were compared to an Objective Performance Goal (OPG) at 12 months.

Results: Procedural success was obtained in 100%.

Results for Doppler ultrasound show primary patency of 87% (77/89, 95% confidence interval 80-94), meeting the requirements of the OPG (73%). 24-month primary patency was 81% (63/78).

The 12-month primary safety endpoint was 96.8% (90/93), meeting the requirements of the OPG (19.6%). 2.2% (2/93) of patients had clinically driven target lesion revascularization, 1.1% (1/93) died, 0 had target limb amputation, 0 had target limb ischemia requiring surgical intervention or surgical repair of the target vessel and 0 had worsening of the Rutherford score by two classes, or class 5 and 6.

The 24-month safety endpoint was 87.8% (79/90). 6.7% (6/90) of patients had clinically driven target lesion revascularization, 3.3% (3/90) died, 0 had target limb amputation, 1.1% (1/90) had target limb ischemia requiring surgical intervention or surgical repair of the target vessel and 1.1% (1/90) had worsening of the Rutherford score by two classes, or class 5 and 6.

14. DIRECTIONS FOR USE

Selection of the device

Refer to the length of the target lesion to establish the length of stent required. In order to appropriately cover the lesion, the deployed stent length should be at least 5 mm longer than the lesion itself on either side.

Refer to the diameter of the vessel proximally and distally of the lesion to establish the appropriate stent size. To ensure a secure placement of the stent, use the larger diameter as the basis for choosing the appropriate stent size. When selecting the model, refer to the information listed in table 1.

Preparing the device

Prior to stent delivery, administer heparin and check that the patient's Activated Clotting Time (ACT) is adequate in conformity with the practice of the Centre.

During preparation of the system, take special care not to let any liquid come into contact with the stent.

- Open the box and take out the bag containing the sterile pouch.
- After examining the package for any damage, open the external non-sterile aluminium foil pouch then remove the content from the internal pouch. (Note: The aluminium foil pouch and the outer surface of the inner plastic pouch **ARE NOT STERILE**. Only the contents of the plastic pouch are STERILE).
- Taking care to place the support on a flat surface in a sterile area, unfasten the band around the handle and remove the catheter from the protective coil.
- Check that the delivery catheter does not present kinks, bends or other damage. Damage to the device may impair its performance.
- Fill a 3 mL syringe with a physiological solution containing heparin, connect to the luer fitting for flushing (figure 1, item 12) and apply positive pressure. Continue flushing until the physiological solution emerges from the distal end of the catheter.
- Inspect the distal end of the catheter, making certain that the stent (figure 1, item 4) is covered completely by the outer sheath.

Do not use the device if damaged or if the stent is partially expanded.

Compatibility with accessories

Introducer sheaths of diameter 6F (2.0 mm) or greater can be used with the NITIDES device.

Guidewires of 0.035 inches (0.89 mm) or less can be used with the NITIDES devices. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

Should it be necessary to utilize a guiding catheter, systems are compatible with guiding catheters of diameter 8F (2.7 mm) or greater. In this instance, make certain that the haemostatic valve is fully open before introducing the selected assembly.

Stent insertion

- Make sure the safety key (figure 1, item 10) of the handle is in the locked position (detail B).
- Wash the exposed portion of the guidewire with a physiological solution containing heparin, to remove traces of blood and contrast medium.
- Advance the system over the guidewire until the stent reaches the treatment site. If unusual resistance is encountered at any time, do not force the system: withdraw the introducer, the guidewire and the device together as a single unit. Applying excessive force and/or incorrect handling may result in the stent being accidentally deployed or in damage to the placement catheter.

Advancing the device

- Advance the system to the point at which the stent is beyond the lesion.
- Draw the system back until the stent is in alignment with the lesion.
- Check that the part of the delivery system outside the body of the patient is freely extended and straight.

Caution: if the delivery system is not properly extended, the stent could be deployed beyond the target lesion.

Stent implantation and expansion

- Firmly grasping the handle with one hand, check that the radiopaque markers on the stent are positioned proximally and distally in relation to the lesion.
- Make certain that the introducer sheath does not move during deployment of the stent.

- Turn the safety key on the handle anticlockwise (figure 1, detail B) to unlock it (STEP 1), then pull it out completely (STEP 2).
- Holding the handle in a fixed position with respect to the patient, check once again that the radiopaque markers are in the desired position with respect to the lesion; begin deploying the stent by turning the wheel (figure 1, item 9) with your thumb in the direction of the arrow, as shown in figure 1.
- Under fluoroscopy, when the distal end of the stent is well-attached to the vessel walls, continue deployment by turning the wheel until the end of the stent is opened and it is completely implanted in the vessel walls.
- If rapid deployment is required, proceed as follows: make sure that there are no obstructions (e.g. female luer fitting screwed) in the rear end of the handle, grasp the rear end of the device (pull-back, figure 1, item 11) with your other hand and pull it, holding the front part of the handle with respect to the patient perfectly still, as shown in figure 1.
- Check fluoroscopically that the stent has fully expanded.
- The stent can be expanded further, if necessary, using a PTA balloon catheter. In selecting the diameter of the balloon catheter and in performing the post-expansion, please refer to the diameter of the vessel that has been evaluated by means of angiography. In any event, take particular care not to expand the stent beyond its nominal diameter, stated on the label.
- The final diameter of the stent must be suitable for the reference vessel diameter.

15. POTENTIAL ADVERSE EVENTS

The possible adverse events that may be associated with the implantation of a peripheral stent include risks related to percutaneous transluminal angioplasty and additional risks associated with the use of stents, as indicated below.

- Access-site infection, pain and haematoma
- Allergic or adverse reaction to contrast medium, drugs used during the procedure or materials composing the stent
- Arterial spasm
- Arterial or stent thrombosis
- Femoral artery aneurysm or pseudoaneurysm
- Atherosclerotic distal embolization
- Bleeding
- Death
- Dissection, perforation or rupture of the stented vessel
- Embolism
- Kidney failure
- Lower limb ischemia (requiring surgical intervention)
- Restenosis of the artery in which the stent was implanted
- Stent embolization
- Vessel occlusion or unresolved stenosis

16. LIABILITY AND WARRANTY

The Manufacturer ensures that the device has been designed, manufactured and packaged with every care using the techniques deemed most suitable among those available at the current state of technology and applying principles of integration of safety in design and construction to ensure safe use when it is used under the conditions and for the intended purposes respecting the precautions indicated in the previous paragraphs, and in any case aimed at reducing as far as possible, but not totally eliminating, the risks related to the use of the device. The device must be exclusively used under the responsibility of specialized medical personnel and taking into account the residual risks and possible side effects and complications of the therapy for which it was intended, including those mentioned in the other sections of this instruction booklet.

The technical complexity and the critical nature of medical choices and methods of application of the device exclude any guarantee of the Manufacturer, explicit or implicit, on the effectiveness of the result following the use of the device, or on its effectiveness in resolving a state of disability, or any other pathology or medical problem. In fact, the results, both from the clinical point of view and in terms of functionality and duration of the device, also depend on numerous factors that are beyond the control of the manufacturer, among which are cited merely as examples and not exhaustive conditions of the patient, the surgical procedure of implantation or use, the methods of manipulation of the device after opening the package.

In light of these factors the manufacturer's liability exists only if it has been ascertained that the device was affected by defects at the time it was put into circulation. For this purpose, the Customer must promptly place the device at the disposal of the Manufacturer, in compliance with the terms of expiry and prescription required by law: the Manufacturer reserves the right to examine the device deemed to be faulty and to establish if it is actually affected by defects that may affect the safety of the product. The manufacturer's liability is excluded if (i) the defect claimed is due to the compliance of the device with a mandatory rule or a binding provision, and / or (ii) the defect that caused the damage did not exist when the Manufacturer released the device and / or (iii) if the state of technical and scientific knowledge at the time the Manufacturer release the device, did not yet allow the device to be considered as defective, and / or (iv) in any other case, also different from those stated above, in which the law expressly excludes, or limits, the responsibility of the Manufacturer for damage from defective products. The guarantee consists exclusively in the replacement of the defective device with another of the same or equivalent type of manufacture of the manufacturer. The warranty applies only if the device is correctly returned to the manufacturer together with a written and detailed report on the claimed defects and, in the case of an implanted device, on the reasons for removal from the patient. The claimed defects must be promptly reported to the Manufacturer in compliance with the terms of forfeiture and statutory prescription.

Once the replacement has been made, the Manufacturer will reimburse the Customer for the expenses incurred by him for the return of the recognized faulty device.

The Manufacturer declines all responsibility in any case of non-observance of the methods of use and of the precautions indicated in this instruction manual and in the case of use after the expiry date indicated on the package.

Furthermore, the Manufacturer declines all responsibility for the consequences of medical choices and the methods of use or application of the device and therefore will not under any circumstances be prosecuted for any material, biological or moral damage of any kind, consequent to the application of the device and to the choice of the implant technique used by the operator.

The agents and representatives of the Manufacturer are not authorized to modify any of the aforementioned conditions nor to make any additional commitment or to make any warranty in relation to this product, other than as indicated above. This warranty is governed by Italian law, unless the mandatory laws of another state are applied, in case the Manufacturer has authorized the circulation of the product in one or more states other than Italy.

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TABLE 1

Usable Catheter Length (UCL)		Unconstrained Stent Diameter (mm)	Vessel Diameter (mm)	Nominal Stent Length (mm)
85 cm	135 cm			
ICNI6020S	ICNI6020L	6.0	4-5	20
ICNI6040S	ICNI6040L	6.0	4-5	40
ICNI6060S	ICNI6060L	6.0	4-5	60
ICNI6080S	ICNI6080L	6.0	4-5	80
ICNI60100S	ICNI60100L	6.0	4-5	100
ICNI60120S	ICNI60120L	6.0	4-5	120
ICNI60150S	ICNI60150L	6.0	4-5	150
ICNI7020S	ICNI7020L	7.0	5-6	20
ICNI7040S	ICNI7040L	7.0	5-6	40
ICNI7060S	ICNI7060L	7.0	5-6	60
ICNI7080S	ICNI7080L	7.0	5-6	80
ICNI70100S	ICNI70100L	7.0	5-6	100
ICNI70120S	ICNI70120L	7.0	5-6	120
ICNI70150S	ICNI70150L	7.0	5-6	150
ICNI8020S	ICNI8020L	8.0	6-7	20
ICNI8040S	ICNI8040L	8.0	6-7	40
ICNI8060S	ICNI8060L	8.0	6-7	60
ICNI8080S	ICNI8080L	8.0	6-7	80
ICNI80100S	ICNI80100L	8.0	6-7	100
ICNI80120S	ICNI80120L	8.0	6-7	120
ICNI80150S	ICNI80150L	8.0	6-7	150