**CiD**  
Carbostent & Implantable Devices

PERIPHERAL VASCULAR

**CFRE**  **BTK**<sup>TM</sup>

AMPHILIMUS<sup>TM</sup> ELUTING PERIPHERAL STENT

Manufactured by: **CiD** S.p.A.  
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy  
**CiD** S.p.A. member of  **Alvimedica** group  
[www.alvimedica.com](http://www.alvimedica.com)

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**CRE8™ BTK**  
**Infrapopliteal Sirolimus-eluting stent coated with „Carbofil™ and rapid-exchange balloon catheter**

**1. DESCRIPTION**

The CRE8™ BTK device consists of a Sirolimus-eluting infrapopliteal stent firmly held on the distal end of a semicompliant balloon catheter.

The stent is a flexible implantable device that can be expanded using a PTA catheter.

The stent is made of cobalt-chromium alloy (L605) and is coated with „Carbofil™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs. Coating the substrate with „Carbofil™ provides it with the bio- and haemocompatible characteristics proper of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

The outer surface of the stent has dedicated grooves, fully coated with „Carbofil™, for containing the pharmaceutical formulation Amphilimus™, which is composed of the drug Sirolimus and a mixture of long-chain fatty acids.

The specific drug dosage for the vessel surface treated is approximately 0.9 µg/mm<sup>2</sup>, which corresponds to a minimum dose of 50 µg on the smaller stent (2.25x8 mm) and a maximum dose of 395 µg on the larger stent (4.0x38 mm-3.5x46 mm).

Two radio-opaque platinum markers at either end of the stent enable accurate positioning on the lesion to be treated.

The rapid-exchange type **balloon catheter** provides a safe means for delivering the infrapopliteal stent to the lesion to be treated.

The distal portion of the catheter has a hydrophilic coating and consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radio-opaque markers, located outside the usable length of the balloon, enable accurate placement across the stenosis.

The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.

Two depth indicators placed at 90 and 100 cm from the distal end, show when the balloon comes out of the guiding catheter.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

The Manufacturer directly manufactures the CRE8™ BTK device and carries out all quality controls both during production and on the finished product, in accordance to the appropriate good manufacturing standards.

**How supplied:**

CRE8™ BTK device is supplied sterile, non-pyrogenic and individually packed in a foil pouch that must not be placed into a sterile field.

**NOTE: The foil pouch is the only sterile barrier.**

The manufacturer sterilization uses an ethylene oxide and CO<sub>2</sub> mixture. Sterility and chemical stability are guaranteed while the package remains intact and correctly stored, until the expiry date printed on the packaging (EXPIRY DATE).

**Contents:**

- One CRE8™ BTK device.
- One compliance table.

**Storage:**

Store in a cool, dry place, away from sunlight, at a temperature of 25°C; excursion permitted to 30 °C.

**2. INTENDED USE**

The stent is indicated to help restore patency to peripheral vessels (e.g. infrapopliteal vessels such as the anterior tibial artery, posterior tibial artery, peroneal artery, tibio-peroneal trunk) in patients suffering from symptomatic chronic ischemia.

**3. INDICATIONS**

The stent is indicated in the following cases:

- Treatment of patients suffering from symptomatic chronic ischemia at infrapopliteal level, characterized by pain of ischemic origin, claudication, ulcers and/or trophic lesions with risk of lower limb amputation, in order to improve or resolve these problems.
- Treatment of obstructive lesions of the infrapopliteal vessels in patients meeting the requirements for a PTA procedure.
- Primary peripheral vessel angioplasty with elective indication of stent implantation in vessels whose diameter is compatible with the models in the catalogue.
- Short- and long-term unsatisfactory results arising from PTA procedures.

**4. CONTRAINDICATIONS**

Use of the stent is contraindicated in the following situations:

- Lesions considered untreatable with PTA or other interventional techniques.
- Disorders / allergies that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to the contrast medium or drugs used during the procedure.
- Lesions on a vessel with a reference diameter < 2.25 mm.
- Patients with known hypersensitivity or allergies to Sirolimus, fatty acids (such as stearic acid, palmitic acid, behenic acid).

**5. WARNINGS**

- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.

- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to the latest literature for information regarding the risks and advantages of the various procedures before choosing the procedure to be adopted.
  - **In patients with known hypersensitivity or allergies to the metal components of the stent, the CRE8™ BTK stent should be used only if the potential benefits of the stent outweigh the potential risks.**
  - There are no adequate or well-controlled studies in pregnant women for this product. The CRE8™ BTK stent should be used during pregnancy only if the potential benefits of the stent outweigh the potential risks.
  - The device is intended for single use only. Do not reuse, reprocess or resterilise. This could compromise device performance and cause contamination risks with the device and cause infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
  - The CRE8™ BTK device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce alterations.
  - The stent should not be touched. This is most important during catheter removal from packaging, placement over the guide wire and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
  - The CRE8™ BTK device is designed to perform as a system. Do not use its components separately.
  - The CRE8™ BTK device is indicated for use in combination with PTA. The delivery catheter is not a peripheral dilatation catheter; it should be used for stent deployment only.
  - Do not use the CRE8™ BTK device if its packaging has been opened or damaged or not correctly stored, or if it is past its expiry date. In such cases product performance and sterility are not guaranteed.
  - Do not use the CRE8™ BTK if the shaft of the catheter has kinks or bends; in such cases, do not try to straighten the catheter.
  - **The stent should not come into contact with liquids prior to preparation and positioning. However, if it is absolutely necessary to irrigate the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).**
  - A surgery team must be available for a possible intervention.
  - Follow the manufacturer's instructions when using accessories (introducer, guiding catheter, guide wire, haemostatic valve).
  - The CRE8™ BTK should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
  - The traditional procedure for stent deployment requires pre-dilation of the lesion.
  - The stent should be implanted in the target lesion using its delivery system.
  - If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the introducer, the guiding catheter and the complete CRE8™ BTK device as a single unit. Applying excessive force and/or incorrectly handling the system may lead to stent deployment or damage to the delivery catheter.
  - If the patient has multiple lesions in a single vessel, it is advisable to first treat the distal lesion, then the proximal lesions. This order will reduce the need to cross the proximal stent during implantation of the distal stent and therefore reduces the risk of displacing the proximal stent.
  - Implantation of an infrapopliteal stent may cause dissection of the vessel distal and/or proximal to the stent and may also cause sudden occlusion of the vessel, making an additional intervention necessary (surgery, further dilatation, placement of additional stents, or other procedures).
  - Do not pressurise the system until the stent is positioned across the lesion to be treated.
  - Never use air or any other gas to inflate the balloon.
  - During catheter inflation do not exceed the Rated Burst Pressure.
  - Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
  - Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
  - If the stent is lost in the vessel, recovery procedures can be initiated. The procedures can, however, cause injury to vessels and/or to the vessel access site.
  - Do not attempt to clean or re-sterilise devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as hazardous medical waste with risk of infection.
  - Procedural success, in terms of restoration of peripheral vascularity, may depend not only on the effectiveness of the device, but also on the characteristics of the treated vessel segment and the downstream perfusion capacity of the district.
  - The physician should inform the patient that if the course of the disease makes the use of bandages and/or physiotherapy necessary, he/she must avoid any local stress to the stent implantation area.
- 6. PRECAUTIONS**
- Before the angioplasty procedure, examine the catheter to check that it operates correctly (no kinks, bends or other damage) and ensure that the dimensions are correct for the specific procedure in which it is to be used.
  - Use of a gauged inflation device is strongly recommended.
  - The CRE8™ BTK device should be used only by physicians trained specifically to perform percutaneous transluminal angioplasty (PTA) and peripheral stent implantation.
  - Where vessels are excessively tortuous and there are also diffuse atheromatous plaques proximal atherosclerosis, the catheter may be difficult to advance. In such cases incorrect handling could lead to dissection or rupture of the vessel.
  - An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter.
  - To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.

## 7. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 2.890 G/cm (28.90 T/m).

### MRI Heating

Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest temperature change	2.4 °C	2.9 °C

### MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

## 8. DRUG REGIME

Studies in clinical literature indicate the need to administer appropriate anticoagulant, anti-platelet and vasodilation therapy for successful stent implantation.

Antiplatelet treatment during the period after the procedure should be used per the guidelines from current ESC or ACC/AHA/SCAI for PCI.

While physicians should adhere to DAPT duration defined into the current ESC (2017 and subsequent update) or ACC/AHA/SCAI Guidelines for PCI, the decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition (e.g. non-deferrable surgery). Guidelines details on indications and modalities for DAPT interruption are available in the ESC or ACC/AHA/SCAI organizations' websites (<https://www.escardio.org/>, <https://www.acc.org/>, <https://www.heart.org/http://www.scai.org/>). In case unanticipated interruption or discontinuation of DAPT is required any time after one month following Cre8™ products family implantation, one-year clinical data of Cre8™ products family indicate low stent thrombosis rates and no observed increased risk for stent thrombosis.

## 9. INTERACTIONS WITH DRUGS

Although no specific clinical data are available, certain drugs, such as Tacrolimus, which act through the same binding protein (FKBP), may interfere with the efficacy of Sirolimus. There have been no studies on the interaction between drugs. Sirolimus is metabolized by CYP3A4. Strong CYP3A4 inhibitors may cause greater exposure to Sirolimus up to levels associated with systemic effects, especially in the case of multiple implants. Systemic exposure to Sirolimus if the patient is taking concomitant systemic immunosuppressive drugs must also be considered.

## 10. POTENTIAL ADVERSE EVENTS

Possible adverse effects of PTA include, but are not limited to:

- Access-site haematoma or pain.
- Arteriovenous fistula.
- Femoral pseudoaneurysm.
- Acute occlusion of the vessel.
- Allergic reaction (to contrast medium, drugs used during the procedure or to the materials composing the stent).
- Death.
- Dissection, perforation, rupture of the artery.
- Distal embolism.
- Infections.
- Restenosis of the vessel.
- Stent occlusion.
- Systemic bleeding.
- Thrombosis (acute, subacute or late).
- Total or partial amputation of the limb.

## 11. DIRECTIONS FOR USE

### Auxiliary material (not included in the box)

- Introducer with a 4F diameter (inner diameter 1.35 mm).
- Guiding catheter(s) with a 5F diameter (inner diameter 1.47 mm) or larger.
- Guidewires 0.014 inches in diameter (0.356 mm) or less.
- Inflation device.
- Syringes.

### Model

Each CRE8™ BTK device is identified by a model code and a lot number; available product codes are listed in Table 1. The code is composed of the letters IC, followed by two letters identifying the type of device, two or three digits indicating the nominal diameter of expansion of the stent and a further two digits indicating the length of the stent.

The batch number allows traceability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels on each box; these labels can be attached to the implant patient's medical chart.

### Stent + delivery catheter preparation

**During preparation of the system, take special care not to let any liquid come into contact with the stent. However, if it is absolutely necessary to irrigate the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).**

The stent should be implanted in the target lesion using its delivery system. Examine the package for any damage and then remove the CRE8™ BTK device and take it into a sterile area.

- Check that the delivery catheter does not have kinks, bends or other damage.
- Carefully slide the protective cover off the stent by grasping the cover at the distal end. Do not touch the part of the cover over or proximal to the stent. Check that the stent is intact and well centred on the balloon.  
**Damage to the CRE8™ BTK device may impair its performance. If the stent is out of place or damaged, do not use it.**
- Flush the guidewire lumen with heparin/saline mixture. **Caution: Do not allow the stent to come into contact with liquids while flushing the guidewire lumen.**
- Prepare the inflation device according to the manufacturer's instructions with dilute contrast medium.
- Purge any air from the balloon the stent is mounted on as follows:
  - Fill the inflation device with 4 ml of diluted contrast medium (50% contrast medium and 50% sterile saline solution).
  - After connecting the inflation device/syringe to the delivery catheter Luer connector, point the distal tip of the catheter (balloon) downwards vertically.
  - Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to neutral as the system fills with contrast medium.**
  - Disconnect the inflation device/ syringe from the Luer connection and remove all air.
  - Without introducing air, repeat steps 2 and 3 and aspirate for 10-15 seconds until air bubbles no longer appear.

### Stent insertion

**Introducers with a 4F diameter (internal diameter 1.35 mm) or larger are compatible for use with the CRE8™ BTK device.**

**Guiding catheters with a 5F diameter (internal diameter 1.47 mm) or larger are compatible for use with the CRE8™ BTK device.**

**Guidewires with a diameter of 0.014 inches (0.356 mm) or less are compatible for use with the CRE8™ BTK device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.**

- Keep the delivery catheter to which the stent is crimped at ambient pressure during insertion (neutral pressure).
- Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- Ensure that the haemostasis valve is completely open before introducing the system.
- Thread the distal tip of the catheter over the proximal end of the guidewire until the guidewire exits the catheter body at the guidewire exit port.
- Slowly advance the dilatation catheter until the appropriate depth marker aligns with the haemostatic valve hub. This indicates that the balloon is exiting the guiding catheter.
- Under fluoroscopy, slowly advance the system over the guidewire to allow the blood retrograde filling of the guiding catheter; keep advancing the system until the stent reaches the treatment site.
- When advancing the delivery system into the target vessel, make sure that the stent and the delivery balloon move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon radiopaque markers in relation to the stent radiopaque markers.

**If resistance is encountered at any time, do not force the system: withdraw the introducer (or the guiding catheter) and the complete CRE8™ BTK device as a single unit. Applying excessive force and/or incorrectly handling the system may lead to stent deployment or damage to the delivery catheter.**

### Stent implantation and expansion

- Position the stent and, with the help of the balloon's radiopaque markers, make sure it is correctly positioned in relation to the lesion to be treated and completely covered.
- Deploy the stent by slowly inflating the balloon to the nominal diameter. Table 2 lists the diameters of the balloon as inflation pressure varies (column I), for the seven nominal diameter groups: 2.25 mm (column II), 2.5 mm (column III), 2.75 mm (column IV), 3.0 mm (column V), 3.5 mm (column VI), 4.0 mm (column VII) and 4.5 mm (column VIII). Values on a dark background refer to pressures above the Rated Burst Pressure\*. The stent expands at a minimum balloon inflation pressure of about 5 atm.
- Stent elastic recoil is between 2 and 7%\* depending on stent model and expansion diameter.**

\* NOTE: These values are obtained from the results of in-vitro testing.

- Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation.**
- Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTA procedure. Ensure the balloon is fully deflated before and remove the catheter. A device with larger and longer balloon requires a longer deflation time.
- If post-dilatation is required, a PTA balloon can be used. Use extreme care: **not to expand the 2.25 mm nominal diameter stent over 2.55 mm. not to expand the 2.50 – 2.75 mm nominal diameter stent over 3.05 mm. not to expand the 3.00 – 3.50 mm nominal diameter stent over 3.85 mm. not to expand the 4.00 – 4.50 mm nominal diameter stent over 5.05 mm.**

**Optimal expansion requires the stent to be in full contact with the arterial wall, so that the stent diameter is equal to the diameter of the reference vessel.**

**ENSURE THAT THE STENT IS FULLY EXPANDED.**

## 12. LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid

conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects.

Under such circumstances, the customer shall deliver the device to the Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty consists exclusively in the replacement of the defective device with another device of the same type or equivalent from the same Manufacturer.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for replacement of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package. Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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SIROLIMUS ELUTING  
CARBOFILM™ COATED PERIPHERAL STENT  
ON RX BALLOON CATHETER


TABLE 1		Stent length							
		8mm	12mm	16mm	20mm	25mm	31mm	38mm	46mm
Nominal stent diameter	2.25 mm	ICLK22508	ICLK22512	ICLK22516	ICLK22520	ICLK22525	ICLK22531	-	-
	2.5 mm	ICLK2508	ICLK2512	ICLK2516	ICLK2520	ICLK2525	ICLK2531	ICLK2538	ICLK2546
	2.75 mm	ICLK27508	ICLK27512	ICLK27516	ICLK27520	ICLK27525	ICLK27531	ICLK27538	ICLK27546
	3.0 mm	ICLK3008	ICLK3012	ICLK3016	ICLK3020	ICLK3025	ICLK3031	ICLK3038	ICLK3046
	3.5mm	ICLK3508	ICLK3512	ICLK3516	ICLK3520	ICLK3525	ICLK3531	ICLK3538	ICLK3546
	4.0mm	ICLK4008	ICLK4012	ICLK4016	ICLK4020	ICLK4025	ICLK4031	ICLK4038	-
	4.5mm	-	ICLK4512	ICLK4516	ICLK4520	ICLK4525	ICLK4531	-	-

TABLE 2		Balloon diameter (mm)							
		2.25	2.50	2.75	3.00	3.50	4.00	4.50	
Inflation pressure (atm)	6	2.08	2.31	2.54	2.77	3.26	3.74	4.15	
	7	2.13	2.37	2.61	2.85	3.35	3.84	4.27	
	8	2.19	2.43	2.68	2.93	3.43	3.92	4.38	
	<b>NP 9</b>	<b>2.25</b>	<b>2.50</b>	<b>2.75</b>	<b>3.00</b>	<b>3.50</b>	<b>4.00</b>	<b>4.50</b>	
	10	2.31	2.55	2.83	3.08	3.58	4.09	4.61	
	11	2.36	2.61	2.90	3.15	3.66	4.17	4.68	
	12	2.40	2.65	2.95	3.20	3.71	4.25	4.74	
	13	2.44	2.70	3.01	3.25	3.77	4.32	4.80	
	14	2.48	2.74	3.06	3.29	3.82	4.37	4.86	
	15	2.51	2.77	3.09	3.32	3.86	4.41	4.92	
	16	2.55	2.80	3.12	3.36	3.90	4.46	4.98	
	17	2.59	2.83	3.15	3.39	3.94	4.50	5.04	
	<b>RBP 18</b>	<b>2.62</b>	<b>2.86</b>	<b>3.19</b>	<b>3.43</b>	<b>3.98</b>	<b>4.54</b>	<b>5.11</b>	
	19	2.65	2.89	3.23	3.47	4.03	4.59	5.19	
	20	2.68	2.93	3.27	3.51	4.08	4.65	5.27	
21	2.71	2.96	3.30	3.55	4.12	4.69	5.33		
22	2.74	2.99	3.34	3.60	4.18	4.74	5.40		

These values result from in vitro testing

- NP Nominal Pressure
- RBP Rated Burst Pressure Do not exceed
- Pressure above Rated Burst Pressure

If postdilatation is required, a PTA balloon can be used. Use extreme care:  
 not to expand the 2.25 mm nominal diameter stents over 2.55 mm  
 not to expand the 2.50 mm and the 2.75 mm nominal diameter stents over 3.05 mm  
 not to expand the 3.00 mm and the 3.50 mm nominal diameter stents over 3.85 mm  
 not to expand the 4.00 mm and the 4.50 mm nominal diameter stents over 5.05 mm