PERIPHERAL VASCULAR

CARBOSTENT™
INPERIA
ADVANCE

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DESCRIPTION

The "INPERIA ADVANCE CARBOSTENT" device consists of an infra-popliteal stent firmly held on the distal end of a semicompliant balloon catheter.

The infra-popliteal **stent** is a flexible implantable device that can be expanded using a PTA balloon catheter.

The stent is made of a cobalt-chromium alloy coated with iCarbofilm™ a thin layer of carbon with a high-density turbostratic structure substantially identical to that of the pyrolytic carbon used for constructing the leaflets of prosthetic mechanical heart valves. Coating the substrate provides it with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Two radio-opaque platinum markers at either end of the stent allow the device to be correctly positioned over the lesion to be treated.

The rapid-exchange type **balloon catheter** provides a safe means for delivering the infra-popliteal stent to the lesion to be treated.

The distal portion of the catheter consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radio-opaque markers, located outside the cylindrical part of the balloon, enable accurate positioning across the stenosis.

The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.

Two depth indicators positioned 90 and 100 cm from the distal end, signal when the balloon exits the guiding catheter for the brachial or femoral approach, respectively.

The proximal end of the catheter has a female Luer lock connector for attachment to the inflation device.

The Manufacturer produces the "INPERIA ADVANCE CARBOSTENT" device directly and carries out all quality controls both during manufacturing and on the finished product, in accordance with the appropriate Good Manufacturing Practices.

INTENDED USE

The stent is intended for use in restoring patency of peripheral vessels (especially infra-popliteal vessels such as the anterior tibial artery, posterior tibial artery, peroneal artery, tibial-peroneal trunk) in patients suffering from symptomatic chronic ischemia.

INDICATIONS FOR USE

The stent is indicated in the following cases:

- Treatment of patients suffering from symptomatic chronic ischemia at an infrapopliteal level, characterized by ischemic pain, claudication, ulcer and/or trophic lesions with the risk of amputation of the lower limbs, in order to improve or heal the above-mentioned problems.
- Treatment of occlusive lesions in infrapopliteal vessels in patients meeting the requirements for a PTA procedure.
- Primary angioplasty of peripheral vessels with elective stenting in vessels having a diameter compatible with the models in the catalogue.
- Unsatisfactory short and long term results originating from PTA procedures

CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations.

- Pregnant women.
- Disorders that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to contrast medium.
- Lesions considered untreatable with PTA or other operating techniques.
- Resistant lesions that cannot be dilated.
- Lesions on a vessel with a reference diameter < 2.25 mm.
- Lesions with severe calcification that resist dilatation with risk of vessel wall rupture.
- Aneurysms immediately adjacent to the stent implant site.

MODEL

Every "INPERIA ADVANCE CARBOSTENT" device is identified by a Model Code and a Lot Number; available models are listed in table 1. The model code is made up of the prefix letters IC followed by two letters indicating the type of the device, two or three digits indicating the nominal expanded diameter of the stent, and two digits indicating the length of the stent.

The lot number allows traceability of all information regarding the manufacturing process and system control in the Manufacturer's Quality Assurance archives.

To make the tracing of the device easier at the user end, the identification data for the stent is printed on the adhesive labels included in each box; these labels can be attached to the patient health record.

PACKAGING

The "INPERIA ADVANCE CARBOSTENT" device is supplied sterile and individually packaged in a single bag that **must not be placed on the sterile field**.

The manufacturer's sterilisation uses an ethylene oxide and carbon dioxide mixture.

Sterility is guaranteed if package is unopened and undamaged and until the use before date printed on the packaging (USE BEFORE).

WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause risks of contamination of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The "INPERIA ADVANCE CARBOSTENT" device is designed as an integral system. Do not use its components separately.
- The "INPERIA ADVANCE CARBOSTENT" device is also indicated for use during endovascular procedures. The delivery catheter is not a peripheral dilatation catheter; it should be used for stent deployment only.
- Do not use the "INPERIA ADVANCE CARBOSTENT" device if its packaging has been opened or damaged or if the period of sterility has expired. In such cases product sterility is not guaranteed.
- Do not remove the stent from the delivery balloon since this could damage the stent and/or cause stent embolization.
- The "INPERIA ADVANCE CARBOSTENT" device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce mechanical alterations.

- Do not use the "INPERIA ADVANCE CARBOSTENT" device if, during handling, due to torsion or resistance to insertion, the proximal portion of the delivery catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The "INPERIA ADVANCE CARBOSTENT" device should only be implanted by physicians specifically trained to perform percutaneous transluminal angioplasty (PTA) and endovascular stent implantation.
- A vascular surgical team must be available for a possible intervention.
- The device should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices. The physician performing the implantation can decide whether to perform the pre-dilatation for the deployment of the stent when there are lesions that permit this approach to be followed.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician can refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- The success of the procedure, in terms of peripheral revascularization, can depend on, along with the effectiveness of the device, the characteristics of the vessel segment being treated and the perfusion capacity of the downstream area.
- The stent should be implanted in the target lesion using its delivery catheter system.
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Do not retract the delivery catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
- Do not attempt to clean or re-sterilise devices that have been in contact with blood or organic tissue. Used devices are to be disposed of as hazardous medical waste with risk of infection.
- Follow the manufacturer's indications for the use of accessories (introducer sheath, guiding catheter, guidewire, haemostatic valve).
- If unusual resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the introducer sheath (or the guiding catheter, if present), the guidewire and the "INPERIA ADVANCE CARBOSTENT" as a single unit. An unexpanded stent can be retracted in the guiding catheter only once. No other movements should be made inside or outside the distal end of the guiding catheter. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.
- If the patient has multiple lesions in one vessel, treat the distal lesion first and then proceed to the proximal lesions. Working in this order reduces the need to cross the proximal stent during implantation of the distal stent, thus lowering the risk of displacing the proximal stent.
- In rare cases, the implanting of a stent may cause dissection of the vessel distal and/or proximal to the implantation site and may also cause sudden occlusion of the vessel, making an additional intervention necessary (surgical vascular replacement operation, further dilatation, placement of additional stents, or other procedures).
- If the stent is lost in the vessel, recovery devices can be used. These procedures can, however, cause injury to vessels and/or to the vascular access site.
- Do not exceed the Rated Burst Pressure during device inflation.
- Never use air or any other gas to inflate the balloon.

PRECAUTIONS

- Be especially careful not to handle the stent on the balloon in any way. This is very important when removing the stent protector, during insertion of the guidewire and when advancing through the haemostatic valve.
- Use of a gauged inflation device is strongly recommended.
- When vessels are excessively tortuous and there is also extensive atheromatous plaque, the catheter may be difficult to advance. In such cases incorrect handling could cause dissection, perforation or rupture of the vessel.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.
- Where multiple stents are needed, the stent materials must be of similar composition.
- For procedural success, it is important that the size of the stent be chosen carefully. As a general rule, the size of the stent must correspond to the diameter of the target vessel and to the length of the lesion. A slightly bigger size is preferred over a smaller one.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
 - Maximum spatial field gradient of 2,890 G/cm (28.90 T/m)
- Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest temperature change	2.4 °C	2.9 °C

In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

DRUG REGIME

Clinical literature indicates the need to administer anticoagulant therapy during the procedure and anti-platelet treatment after the procedure.

Given the complexity of the procedure and the many factors influencing patient condition, definition of the appropriate therapy is left to the discretion of the physician.

INSTRUCTIONS FOR USE

Stent + delivery catheter system preparation

Before commencing the angioplasty procedure, administer heparin and check that the patient's Activated Clotting Time (ACT) is greater than 300 seconds.

The stent should be implanted in the target lesion using the delivery catheter system.

- After having examined the package for any damage, remove the "INPERIA ADVANCE CARBOSTENT" device in a sterile area.
- Check that the delivery catheter does not present kinks, bends or other damage.
- Carefully slide the protector off the stent. Check that the stent is not damaged and is in the centre of the balloon.

Damage to the "INPERIA ADVANCE CARBOSTENT" device may impair its performance. Do not use the stent if it is displaced or damaged.

- Flush the guidewire lumen with heparinized saline solution.
- Prepare the inflation device according to the manufacturer's instructions.
- Purge any air from the balloon on which the stent is mounted as follows:
 - Fill the inflation device with approximately 4 ml of contrast medium;
 - After connecting the inflation device to the delivery catheter's Luer connector, point the distal tip of the balloon downwards.
 - Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to normal**, as the system fills with contrast medium.
- Without introducing air, repeat this operation and aspirate for 10-15 seconds until bubbles no longer appear.

Stent insertion

Introducer sheaths of diameter 4F (ID 1.35 mm) or greater can be used with the "INPERIA ADVANCE CARBOSTENT" devices.

Guiding catheters of diameter 5F (ID 1.47 mm) or greater can be used with the "INPERIA ADVANCE CARBOSTENT" devices.

Guidewires having a diameter of 0.014 inches (0.356 mm) or less can be used with the "INPERIA ADVANCE CARBOSTENT" device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- Keep the delivery catheter on which the stent is pre-mounted at ambient pressure during insertion.
- Wash the exposed portion of the guidewire with a physiological solution containing heparin in order to remove traces of blood and contrast medium.
- Ensure that the haemostatic valve is completely open before introducing the stent + delivery catheter system.
- Slowly advance the system over the guidewire to allow the retrograde filling of the guiding catheter with blood; continue advancing the system until the stent reaches the target site.
- While advancing the system, fluoroscopy should be used to ensure that the stent and balloon move as a single unit by checking the position of the respective radio-opaque markers.

CAUTION: If unusual resistance is encountered at any time, do not force the system: withdraw the introducer sheath (or guiding catheter) and the stent + delivery catheter system as if it were a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.

Implanting and expansion of the stent

- Position the stent across the target lesion.
- Expand the stent by slowly inflating the balloon to the nominal diameter. Table 2 lists the diameters of the system as balloon inflation pressure varies (column I), for the seven nominal diameter groups: 2.25 mm (column II), 2.5 mm (column III), 2.75 mm (column IV), 3.0 mm (column V), 3.5 mm (column VI), 4.0 mm (column VII), 4.5 mm (column VIII). Values on a dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a minimum balloon inflation pressure of about 6 atm.

The stent elastic recoil is between 2 and 7%* depending on stent model and expansion diameter.

* NOTE: These values were obtained through in-vitro testing

- Deflate the balloon and angiographically check that the stent has fully expanded. If necessary, re-inflate the balloon in order to obtain the optimal diameter.
- Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTA procedure.
- If post-dilatation is necessary, a PTA balloon can be used. Use extreme care:
 - not to expand the 2.25 mm nominal diameter stent past 2.55 mm.
 - not to expand the 2.50 mm nominal diameter stent past 3.05 mm.
 - not to expand the 2.75 mm nominal diameter stent past 3.05 mm.
 - not to expand the 3.00 mm nominal diameter stent past 3.85 mm.
 - not to expand the 3.50 mm nominal diameter stent past 3.85 mm.
 - not to expand the 4.00 mm nominal diameter stent past 5.05 mm.
 - not to expand the 4.50 mm nominal diameter stent past 5.05 mm.

The final diameter of the stent must be suitable for the diameter of the reference vessel.

Ensure that the stent is fully expanded.

COMPLICATIONS

Stent implantation may give rise to the following complications:

- Death
- Dissection, perforation, rupture and/or damage of the infrapopliteal artery
- Distal and/or stent embolism
- Acute vessel occlusion (acute or sub-acute thrombosis localised at stent level)
- Chronic vessel occlusion (restenosis of the lesion at stent level)
- Arterial spasm
- Infection
- Access site haematoma
- Access site bleeding
- Drug-induced side effects (anti-coagulant/anti-platelet)
- Contrast medium-induced side effects
- Acute or sub-acute thrombosis localised at stent level
- Total or partial limb amputation
- Artero-venous fistula
- Femoral pseudoaneurysm

- MRI artifacts
- Systemic haemorrhage
- Allergy development

LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the utmost care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforementioned conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. Under such circumstances, the customer shall deliver the device to The Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty exclusively involves the replacement of the defective device with another device of the same type or equivalent produced by the Manufacturer.

The warranty only applies if the device is returned correctly packaged to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for replacement of the device proven to be defective.

The Manufacturer declines any responsibility for cases of failing to observe the methods of use and the precautions stated in this instructions booklet and for cases where the device is used after the use before printed on the package.

Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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INPERIA

ADVANCE

Carbofilm™ coated peripheral infrapopliteal CoCr Stent on RX Balloon Catheter

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Pressure above Rated Burst Pressure
Do not exceed Rated Burst Pressure

These values result from in vitro testing

If postdilatation is required, a PTA balloon can be used. Use extreme care:
not to expand the 2.25 mm nominal diameter stents over 2.55 mm
not to expand the 2.50 mm and the 2.75 mm nominal diameter stents over 3.05 mm
not to expand the 3.00 mm and the 3.50 mm nominal diameter stents over 3.85 mm
not to expand the 4.00 mm and the 4.50 mm nominal diameter stents over 5.05 mm